



## SEAMLESS COMMUNICATION STARTS WITH SOUNDGATE

Convenience. Connectivity. Comfort.



[www.bernafon.com](http://www.bernafon.com)


**bernafon**<sup>®</sup>  
Your hearing · Our passion

## Go Wireless, Relish the Convenience

You are happy with your hearing instruments. They cater to your outgoing nature and active lifestyle as you go to the movies, enjoy live music, and engage in table conversation.

Voices are clear and natural sounding. Optimal speech understanding and listening comfort have significantly improved your quality of life.

And yet, in our communication-technology-driven times, you sometimes wish your hearing instruments also worked like ...



### ... Like a Wireless Headset

This is exactly what SoundGate does. Our easy-to-use interface between hearing instruments and Bluetooth®-enabled cellular phones lets you:

- talk wirelessly and hands-free – enjoy greater mobility;
- listen to music via cell phone or MP3 player – enjoy greater flexibility;
- manage volume and program controls – enjoy greater discretion.


**CONVENIENCE**

## One Device. Countless Options. At the Push of a Button.

Pair SoundGate with all your external audio sources. It is compatible with all wireless-enabled, behind-the-ear Bernafon hearing systems.



**Comfort**  
The TV adapter enables wireless reception of TV audio signals.



**Convenience**  
The phone adapter provides hands-free landline phone calls.



**Mobility**  
Pair SoundGate with Bluetooth®-enabled audio devices – including cell phones, computers, and MP3 players – and enjoy uninterrupted audio streaming in both ears simultaneously. At home, at work, and on the go.



# CONNECTIVITY

## Shared TV Enjoyment

Small and portable, the TV adapter easily connects to any television set.

Since TV signals are directly transmitted to your hearing instruments via SoundGate, you can watch your favorite shows at your preferred volume level. And so can your friends and family.

With its AutoConnect feature, the TV adapter is always ready for use within 10 meters (30 feet) of SoundGate.

Worried about missing a call? Don't. SoundGate alerts you of incoming calls on your cell or ...



## ... Your Landline Phone

SoundGate not only enables hands- and hassle-free calls from your cellular, but also from your landline phone.

Paired with our phone adapter, SoundGate transmits the signals from the stationary phone straight to your hearing instruments via Bluetooth® wireless technology.

Automatically ready for use within 10 meters (30 feet) of SoundGate, the phone adapter delivers uninterrupted landline phone calls in high fidelity.

**To learn more about SoundGate contact your local hearing care professional today.**

COMFORT

Hearing and understanding are basic human needs. Hearing problems and the resulting social consequences represent a challenge to society. We believe in a world where people with restricted hearing can communicate again thanks to advanced technology.

For more information about hearing and hearing loss [www.bernafon.com](http://www.bernafon.com)

**SWISS**   
Engineering

**bernafon**   
*Your hearing · Our passion*

**Berna fon Companies:** Australia · Canada · Denmark · Finland · France · Germany · Italy · Japan · Korea · Netherlands · New Zealand · Poland · Sweden · Switzerland · UK · USA